DiagnosticaVET

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Anaplasma and ehrlichia are vector-borne (ticks), intracellular bacteria. *Dirofilaria immitis* and the **dog's leishmaniasis** are transmitted by blood-sucking mosquitoes.

D. immitis (CHW): Infection with *D. immitis* can lead to a serious and even fatal disease in dogs. Adult heartworms can survive in dogs for 5-7 years. The antigen test is used to detect circulating *D. immitis* antigens in the blood; adult female heartworms are detected.

Anaplasmosis: In Europe, **Anaplasma phagocytophilum** is found in domestic dogs. The antibodies are detected using *A. phagocytophilum* antigens.

Ehrlichiosis: Dogs infected with *Ehrlichia* develop canine monocytic ehrlichiosis. The occurrence of infection is related to the spread of the vector *R. sanguineus*. The antibodies are detected using *E. canis* antigens.

Leishmaniasis: In Europe, **dog leishmaniasis** is caused by the protozoon *Leishmania infantum*. The antibodies are detected by a specific antigen and can be used in diseased dogs about six to eight weeks after initial infection.

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
Remarks for the test result interpretation:

Negative test results in **antibody tests**: they do not rule out an infection, as seronegative results can occur at any stage of infection.

Positive test result for antibody tests: In the endemic area, positive antibody results can be due to a previous infection.

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SIMBOLI

	Only for one use		Read user instruction carefully
	Content		Storage temperature
	Lot number		Expiry date

Disclaimer

The Entire Risk Of Liability In Connection With The Use Of This Product Is With The Buyer. The Manufacturer Is Not Liable For Any Indirect, Special Or Consequential Damages Of Any Sort Resulting From Use Of This Product And The Test Evaluation

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